


(S)-Y-20811

Chemical Properties

CAS No. : 122331-75-3
 Formula: C₂₀H₁₉N₂NaO₃
 Molecular Weight: 358.37
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Y-20811, (S)-, is a compound that inhibits thromboxane A ₂ synthase.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7904 mL	13.9521 mL	27.9041 mL
5 mM	0.5581 mL	2.7904 mL	5.5808 mL
10 mM	0.279 mL	1.3952 mL	2.7904 mL
50 mM	0.0558 mL	0.279 mL	0.5581 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481