


Truxicurium iodide

Chemical Properties

CAS No. : 4304-01-2
 Formula: C₃₄H₅₂I₂N₂O₄
 Molecular Weight: 806.6
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Truxicurium iodide is a neuromuscular blocking agent.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.2398 mL	6.1989 mL	12.3977 mL
5 mM	0.248 mL	1.2398 mL	2.4795 mL
10 mM	0.124 mL	0.6199 mL	1.2398 mL
50 mM	0.0248 mL	0.124 mL	0.248 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481