


WAY-313013

Chemical Properties

CAS No. : 326898-76-4
 Formula: C₁₂H₁₄BrNO₂
 Molecular Weight: 284.15
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	WAY-313013 is an inhibitor specifically targeting p38a Mitogen-Activated Protein Kinase (MAPK).
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.5193 mL	17.5963 mL	35.1927 mL
5 mM	0.7039 mL	3.5193 mL	7.0385 mL
10 mM	0.3519 mL	1.7596 mL	3.5193 mL
50 mM	0.0704 mL	0.3519 mL	0.7039 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481