


## Biclodil

## Chemical Properties

CAS No. :	85125-49-1
Formula:	C <sub>8</sub> H <sub>8</sub> Cl <sub>2</sub> N <sub>4</sub> O
Molecular Weight:	247.08
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Research indicates that Biclodil is a vasodilator with blood pressure-lowering effects.
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0473 mL	20.2364 mL	40.4727 mL
5 mM	0.8095 mL	4.0473 mL	8.0945 mL
10 mM	0.4047 mL	2.0236 mL	4.0473 mL
50 mM	0.0809 mL	0.4047 mL	0.8095 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481