

CB-03-10

Chemical Properties


CAS No. : 1902166-22-6

Formula: C₂₉H₄₂O₆

Molecular Weight: 486.64

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CB-03-10 is a dual inhibitor of AR/GR.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0549 mL	10.2745 mL	20.5491 mL
5 mM	0.411 mL	2.0549 mL	4.1098 mL
10 mM	0.2055 mL	1.0275 mL	2.0549 mL
50 mM	0.0411 mL	0.2055 mL	0.411 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481