


Spiclamine

Chemical Properties

CAS No. :	90243-97-3
Formula:	C ₂₀ H ₂₅ ClN ₂ O
Molecular Weight:	344.88
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Spiclamine (Epsilonamine) is a hemostatic agent.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8996 mL	14.4978 mL	28.9956 mL
5 mM	0.5799 mL	2.8996 mL	5.7991 mL
10 mM	0.290 mL	1.4498 mL	2.8996 mL
50 mM	0.058 mL	0.290 mL	0.5799 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481