


SM-8849

Chemical Properties

CAS No. : 113759-19-6
 Formula: C₁₈H₁₇FN₂S
 Molecular Weight: 312.4
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	SM-8849 is a novel thiazole derivative with anti-arthritis activity.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.201 mL	16.0051 mL	32.0102 mL
5 mM	0.6402 mL	3.201 mL	6.402 mL
10 mM	0.3201 mL	1.6005 mL	3.201 mL
50 mM	0.064 mL	0.3201 mL	0.6402 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481