

TA 078

Chemical Properties


CAS No. : 85915-86-2

Formula: C₁₈H₂₃NO₄ 1/2C₄H₄O₄

Molecular Weight: 375.415

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	TA 078 is an antidiabetic compound that effectively suppresses the rise in blood glucose levels when administered orally to mice undergoing glucose load.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6637 mL	13.3184 mL	26.6368 mL
5 mM	0.5327 mL	2.6637 mL	5.3274 mL
10 mM	0.2664 mL	1.3318 mL	2.6637 mL
50 mM	0.0533 mL	0.2664 mL	0.5327 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481