


Furanochromone

Chemical Properties

CAS No. : 67777-61-1
Formula: C₁₁H₆O₃
Molecular Weight: 186.16
Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	FURANOCHROMONE is a type of plant-based furan compound.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.3717 mL	26.8586 mL	53.7172 mL
5 mM	1.0743 mL	5.3717 mL	10.7434 mL
10 mM	0.5372 mL	2.6859 mL	5.3717 mL
50 mM	0.1074 mL	0.5372 mL	1.0743 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481