

LR132

Chemical Properties


CAS No. : 135211-15-3

Formula: C₁₈H₂₆Cl₂N₂

Molecular Weight: 341.32

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

LR132 is a novel sigma receptor ligand with antagonistic activity.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9298 mL	14.649 mL	29.298 mL
5 mM	0.586 mL	2.9298 mL	5.8596 mL
10 mM	0.293 mL	1.4649 mL	2.9298 mL
50 mM	0.0586 mL	0.293 mL	0.586 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481