


Texacromil

Chemical Properties

CAS No. : 77005-28-8
 Formula: C₁₄H₁₄O₆S
 Molecular Weight: 310.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Texacromil, patented by C. M. Industries S. A., is a benzochromene derivative used as an inhibitor of passive cutaneous anaphylaxis.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2225 mL	16.1124 mL	32.2248 mL
5 mM	0.6445 mL	3.2225 mL	6.445 mL
10 mM	0.3222 mL	1.6112 mL	3.2225 mL
50 mM	0.0644 mL	0.3222 mL	0.6445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481