


RG-7795

Chemical Properties

CAS No. : 892498-64-5
 Formula: C₁₂H₁₄N₄O₅S
 Molecular Weight: 326.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	RG-7795 is a small molecule prodrug whose active component functions as a toll-like receptor 7 (TLR 7) agonist.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0644 mL	15.3219 mL	30.6438 mL
5 mM	0.6129 mL	3.0644 mL	6.1288 mL
10 mM	0.3064 mL	1.5322 mL	3.0644 mL
50 mM	0.0613 mL	0.3064 mL	0.6129 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481