

## Phenampromide tartrate

## Chemical Properties


CAS No. : 7645-04-7

Formula: C<sub>17</sub>H<sub>26</sub>N<sub>2</sub>O x C<sub>4</sub>H<sub>6</sub>O<sub>6</sub>

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description

Phenampromide tartrate is an opioid analgesic with a structure similar to that of isomethadone.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E\_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481