

VUF11418

## Chemical Properties


CAS No. : 1414376-85-4

Formula: C<sub>25</sub>H<sub>31</sub>I<sub>2</sub>N

Molecular Weight: 599.33

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	VUF11418 is a selective non-peptide mimetic agonist targeting CXCR3.
Targets(IC50)	CXCR

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6685 mL	8.3426 mL	16.6853 mL
5 mM	0.3337 mL	1.6685 mL	3.3371 mL
10 mM	0.1669 mL	0.8343 mL	1.6685 mL
50 mM	0.0334 mL	0.1669 mL	0.3337 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481