


Triclonide

Chemical Properties

CAS No. : 26849-57-0
 Formula: C₂₄H₂₈Cl₃FO₄
 Molecular Weight: 505.83
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Triclonide is an anti-inflammatory steroid for topical use.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.9769 mL	9.8847 mL	19.7695 mL
5 mM	0.3954 mL	1.9769 mL	3.9539 mL
10 mM	0.1977 mL	0.9885 mL	1.9769 mL
50 mM	0.0395 mL	0.1977 mL	0.3954 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481