


IRAK4 ligand-12

Chemical Properties

CAS No. :	2573304-97-7
Formula:	C23H24F2N4O3
Molecular Weight:	442.46
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	IRAK4 ligand-12 is a target protein ligand for PROTAC (Ligand for Target Protein for PROTAC) and is utilized in the synthesis of KTX-951.
Targets(IC50)	Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2601 mL	11.3005 mL	22.6009 mL
5 mM	0.452 mL	2.2601 mL	4.5202 mL
10 mM	0.226 mL	1.130 mL	2.2601 mL
50 mM	0.0452 mL	0.226 mL	0.452 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481