

hu14.18 (K322A)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	hu14.18 (K322A) is a humanized monoclonal antibody targeting disialoganglioside GD2, comprised of a fully human amino acid sequence of immunoglobulin G1 heavy chain and κ light chain. This antibody exhibits antibody-dependent cellular cytotoxicity (ADCC) and is utilized in studying high-risk neuroblastoma.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481