


## 1,2-Dielaidoyl-3-palmitoyl-rac-glycerol

## Chemical Properties

CAS No. :	68069-77-2
Formula:	C55H102O6
Molecular Weight:	859.4
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	1,2-Dielaidoyl-3-palmitoyl-rac-glycerol is a triglyceride featuring Elaidic acid at both the sn-1 and sn-2 positions, and Palmitic acid at the sn-3 position.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1636 mL	5.818 mL	11.636 mL
5 mM	0.2327 mL	1.1636 mL	2.3272 mL
10 mM	0.1164 mL	0.5818 mL	1.1636 mL
50 mM	0.0233 mL	0.1164 mL	0.2327 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481