


Sporidesmolide II

Chemical Properties

CAS No. :	3200-75-7
Formula:	C ₃₄ H ₆₀ N ₄ O ₈
Molecular Weight:	652.86
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Sporidesmolide II is a cyclic depsipeptide that was initially isolated from <i>P. chartarum</i> .
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.5317 mL	7.6586 mL	15.3172 mL
5 mM	0.3063 mL	1.5317 mL	3.0634 mL
10 mM	0.1532 mL	0.7659 mL	1.5317 mL
50 mM	0.0306 mL	0.1532 mL	0.3063 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481