


PTX80

Chemical Properties

CAS No. : 2376297-69-5
 Formula: C₂₆H₂₆N₄O₃S
 Molecular Weight: 474.58
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	PTX80 is an antagonist of p62, with an IC ₅₀ value of 31.18 nM. It has been shown to reduce tumor volume in a HCT116 colorectal cancer mouse xenograft model.
Targets(IC ₅₀)	Others,p62

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.1071 mL	10.5356 mL	21.0713 mL
5 mM	0.4214 mL	2.1071 mL	4.2143 mL
10 mM	0.2107 mL	1.0536 mL	2.1071 mL
50 mM	0.0421 mL	0.2107 mL	0.4214 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481