


Torososide B

Chemical Properties

CAS No. :	243835-63-4
Formula:	C40H52O25
Molecular Weight:	932.83
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Torososide B is an anti-allergic compound that inhibits leukotriene release from calcium ionophore A-stimulated rat peritoneal mast cells.
Targets(IC50)	Leukotriene Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.072 mL	5.360 mL	10.7201 mL
5 mM	0.2144 mL	1.072 mL	2.144 mL
10 mM	0.1072 mL	0.536 mL	1.072 mL
50 mM	0.0214 mL	0.1072 mL	0.2144 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481