

Dextran-CY5.5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran-CY5.5 is a conjugate of Dextran labeled with CY5.5. CY5.5 is a cyanine dye (with excitation and emission wavelengths typically at 673 nm and 707 nm, respectively) utilized for labeling proteins, antibodies, and small molecule compounds. Dextran inhibits platelet aggregation and coagulation factors, serving as a plasma volume expander.
Targets(IC50)	Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481