

DEAE-dex-Cy5

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	DEAE-dex-Cy5 is a fluorescent labeling reagent that combines the Cy5 fluorescent dye with Dextran (dex) and diethylaminoethyl (DEAE). Cy5, a cyanine dye, has a wavelength of approximately 662 nm and is used for molecular labeling and microscopic imaging. Dextran inhibits platelet aggregation and coagulation factors, serving as a plasma volume expander.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481