


Pyrenophorin

Chemical Properties

CAS No. : 5739-85-5
 Formula: C₁₆H₂₀O₆
 Molecular Weight: 308.33
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	Pyrenophorin ((-)-Pyrenophorin) is a fungal metabolite found in <i>D. avenae</i> with phytotoxic properties. It exhibits potent cytotoxicity against several cancer cell lines, with IC ₅₀ values ranging from 0.07 to 7.8 μM.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2433 mL	16.2164 mL	32.4328 mL
5 mM	0.6487 mL	3.2433 mL	6.4866 mL
10 mM	0.3243 mL	1.6216 mL	3.2433 mL
50 mM	0.0649 mL	0.3243 mL	0.6487 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481