

Dextran-CY3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Dextran-CY3 is a fluorescent labeling reagent combining the Cy3 fluorescent dye with Dextran. The Cy3 fluorophore is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. Cy3 has a maximum emission wavelength of approximately 562-570 nm. Dextran inhibits platelet aggregation and clotting factors, and is used as a plasma volume expander.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481