

E3 Ligase Ligand-linker Conjugate 141

Chemical Properties

CAS No. : 2229724-20-1


Formula: C₂₀H₂₄N₄O₃

Molecular Weight: 368.43

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	E3 Ligase Ligand-linker Conjugate 141 is a combination of an E3 ubiquitin ligase ligand and a linker, utilized in the synthesis of PROTACER Degradator-12.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7142 mL	13.5711 mL	27.1422 mL
5 mM	0.5428 mL	2.7142 mL	5.4284 mL
10 mM	0.2714 mL	1.3571 mL	2.7142 mL
50 mM	0.0543 mL	0.2714 mL	0.5428 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481