

SGN-B6A

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
-20°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	SGN-B6A is an ADC (antibody-drug conjugate) that targets integrin $\beta$ -6 (ITGB6) using the human IgG1 monoclonal antibody Sigvotatug, and exhibits cytotoxicity towards various integrin $\beta$ -6 positive cancer cells through the mitosis inhibitor MMAE.
Targets(IC50)	Antibody-Drug Conjugates (ADCs),Integrin

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481