


Glucoberteroin

Chemical Properties

CAS No. :	29611-01-6
Formula:	C ₁₃ H ₂₅ NO ₉ S ₃
Molecular Weight:	435.53
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Glucoberteroin is a type of thioglycoside that can be isolated from Degenia velebitica. It inhibits protein glycation and carbonylation, thereby reducing the formation of carbonyl proteins, which aids in preventing skin aging.
Targets(IC50)	Endogenous Metabolite

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2961 mL	11.4803 mL	22.9605 mL
5 mM	0.4592 mL	2.2961 mL	4.5921 mL
10 mM	0.2296 mL	1.148 mL	2.2961 mL
50 mM	0.0459 mL	0.2296 mL	0.4592 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481