

## Alexa Fluor 430 NHS ester

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alexa Fluor 430 NHS ester is a compound used to conjugate the green fluorescent dye Alexa Fluor 430 with proteins or antibodies. Alexa Fluor 430 provides stable signal generation in imaging and flow cytometry. The NHS ester is suitable for labeling proteins, amine-modified oligonucleotides, and other primary amine-containing molecules (R-NH <sub>2</sub> ).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481