

Anti-DM1 antibody(14E3) (Rabbit mAb)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
-20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-DM1 antibody (14E3) (Rabbit mAb) is a rabbit-derived IgG antibody targeting DM1, can be used for activity assays, functional studies, binding affinity analysis, and other applications, serving as an essential tool in biological research.
Targets(IC50)	Microtubule Associated,ADC Antibody

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481