


4-Bromo-2,5-DMMA

Chemical Properties

CAS No. :	155638-80-5
Formula:	C ₁₂ H ₁₈ BrNO ₂
Molecular Weight:	288.18
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	4-Bromo-2,5-DMMA (Compound 2) demonstrates affinity for the 5-HT ₂ binding site and has an ED ₅₀ of 0.82 mg/kg in rat discrimination experiments.
Targets(IC ₅₀)	5-HT Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4701 mL	17.3503 mL	34.7005 mL
5 mM	0.694 mL	3.4701 mL	6.9401 mL
10 mM	0.347 mL	1.735 mL	3.4701 mL
50 mM	0.0694 mL	0.347 mL	0.694 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481