


(S)-3-Aminopyrrolidine-2,5-dione

Chemical Properties

CAS No. :	73537-92-5
Formula:	C ₄ H ₆ N ₂ O ₂
Molecular Weight:	114.1
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(S)-3-Aminopyrrolidine-2,5-dione is a compound with potential anticonvulsant properties, and its derivatives exhibit significant in vivo efficacy.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	8.7642 mL	43.8212 mL	87.6424 mL
5 mM	1.7528 mL	8.7642 mL	17.5285 mL
10 mM	0.8764 mL	4.3821 mL	8.7642 mL
50 mM	0.1753 mL	0.8764 mL	1.7528 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481