


BDP TMR carboxylic acid

Chemical Properties

CAS No. :	287384-28-5
Formula:	C ₂₁ H ₂₁ BF ₂ N ₂ O ₃
Molecular Weight:	398.21
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <i>Actual storage temperature shall be subject to the COA.</i>



Biological Description

Description	BDP TMR carboxylic acid is a BDP TMR-labeled carboxylic acid, featuring the BDP fluorophore. It is suitable for labeling oligonucleotides and sequencing amino acids.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5112 mL	12.5562 mL	25.1124 mL
5 mM	0.5022 mL	2.5112 mL	5.0225 mL
10 mM	0.2511 mL	1.2556 mL	2.5112 mL
50 mM	0.0502 mL	0.2511 mL	0.5022 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481