


## Calcium Channel antagonist 5

### Chemical Properties

CAS No. : 21829-09-4  
 Formula: C<sub>17</sub>H<sub>18</sub>N<sub>2</sub>O<sub>6</sub>  
 Molecular Weight: 346.34  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*



### Biological Description

Description	Calcium Channel antagonist5 (compound 32) is a calcium channel antagonist with a pIC <sub>50</sub> of 5.50.
Targets(IC <sub>50</sub> )	Calcium Channel

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.8873 mL	14.4367 mL	28.8734 mL
5 mM	0.5775 mL	2.8873 mL	5.7747 mL
10 mM	0.2887 mL	1.4437 mL	2.8873 mL
50 mM	0.0577 mL	0.2887 mL	0.5775 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481