


20-Hydroxyganoderic acid A

Chemical Properties

CAS No. :	2374206-89-8
Formula:	C ₃₀ H ₄₄ O ₈
Molecular Weight:	532.67
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	20-Hydroxyganoderic acid A is a triterpene compound identified in Ganoderma lucidum. This compound may exhibit anti-inflammatory properties, suggesting its potential application in the research of inflammatory diseases.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8773 mL	9.3867 mL	18.7733 mL
5 mM	0.3755 mL	1.8773 mL	3.7547 mL
10 mM	0.1877 mL	0.9387 mL	1.8773 mL
50 mM	0.0375 mL	0.1877 mL	0.3755 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481