

5-FAM-LL-37 TFA

## Chemical Properties

CAS No. :

Formula: C226H350N60O59.xC2HF3O2

Molecular Weight:

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	5-FAM-LL-37 TFA is the TFA salt form of 5-FAM-LL-37. This fluorescein-labeled LL-37 peptide maintains the antimicrobial and immunomodulatory properties of LL-37. It binds to bacterial cell membranes, disrupting their integrity, and exhibits broad-spectrum antibacterial activity.
Targets(IC50)	Antibacterial

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481