

P-ESBP-DOX

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	P-ESBP-DOX is an HPMA copolymer-drug conjugate composed of an E-selectin binding peptide and Doxorubicin. It exhibits cytotoxicity against TNF α -activated human vascular endothelial cells (IVEC) with an IC50 of 0.28 μ M and is useful for tumor vasculature research.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481