


KCC2 potentiators-1

Chemical Properties

CAS No. :	3058906-46-7
Formula:	C ₁₅ H ₁₇ N ₃ O ₅
Molecular Weight:	287.38
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Compound 2, known as KCC2 potentiators-1, acts as an enhancer of the potassium chloride cotransporter-2 (KCC2). It holds potential for applications in studies related to neurological disorders.
Targets(IC50)	Others,Na-K-Cl cotransporter

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.4797 mL	17.3986 mL	34.7971 mL
5 mM	0.6959 mL	3.4797 mL	6.9594 mL
10 mM	0.348 mL	1.7399 mL	3.4797 mL
50 mM	0.0696 mL	0.348 mL	0.6959 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481