

## DOX-PEG-Cy3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	DOX-PEG-Cy3 (Doxorubicin-PEG-Cy3) is a Cy3-labeled DOX-PEG conjugate. The Cy3 fluorophore is commonly used for applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry. The maximum emission wavelength of Cy3 is approximately 562-570 nm. DOX is a broad-spectrum anthracycline antibiotic with cytotoxic properties.
Targets(IC50)	Antibiotic

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481