


Lenalidomide-Glycolic acid

Chemical Properties

CAS No. :	2229976-15-0
Formula:	C ₁₉ H ₂₂ N ₂ O ₆
Molecular Weight:	374.39
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Lenalidomide-Glycolic acid is a conjugate of an E3 ligase ligand and a linker, used in the synthesis of the PROTAC-type FLT3 degrader LWY713.
Targets(IC50)	E3 Ligase Ligand-Linker Conjugates

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.671 mL	13.3551 mL	26.7101 mL
5 mM	0.5342 mL	2.671 mL	5.342 mL
10 mM	0.2671 mL	1.3355 mL	2.671 mL
50 mM	0.0534 mL	0.2671 mL	0.5342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481