


## APP-FUBINACA

## Chemical Properties

CAS No. :	1185282-03-4
Formula:	C <sub>24</sub> H <sub>21</sub> FN <sub>4</sub> O <sub>2</sub>
Molecular Weight:	416.45
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	APP-FUBINACA is a cannabinoid receptor agonist, classified as a phenylalanine amide-based indazole-3-carboxamide derivative. It exhibits neurostimulatory effects.
Targets(IC50)	Cannabinoid Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4012 mL	12.0062 mL	24.0125 mL
5 mM	0.4802 mL	2.4012 mL	4.8025 mL
10 mM	0.2401 mL	1.2006 mL	2.4012 mL
50 mM	0.048 mL	0.2401 mL	0.4802 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481