

## Dextran-CY5

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Dextran-CY5 is dextran labeled with the fluorescent dye CY5, which has a peak excitation wavelength of 633 nm and a peak emission wavelength of 670 nm. Dextran-CY5 is useful for visualizing the specific subcellular distribution of glucans at the nanoscale level.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481