

L321

Chemical Properties

CAS No. :

Formula: C₂₃H₂₂F₃N₄O₄

Molecular Weight: 433.42

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	4-Methylphenoxyacetic acid-Oxindole-CF ₃ is the E3 ligase ligand component of L134. It is suitable for use in PROTAC synthesis.
Targets(IC ₅₀)	Ligands for E3 Ligase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3072 mL	11.5362 mL	23.0723 mL
5 mM	0.4614 mL	2.3072 mL	4.6145 mL
10 mM	0.2307 mL	1.1536 mL	2.3072 mL
50 mM	0.0461 mL	0.2307 mL	0.4614 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481