

ITK ligand 1

Chemical Properties

CAS No. :

Formula: C₂₂H₂₆ClN₅O

Molecular Weight: 411.93

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ITK ligand 1 is a PROTAC target protein ligand (Ligands for Target Protein for PROTACS) that is used in the synthesis of PROTACS.
Targets(IC50)	Ligands for Target Protein for PROTAC, Tyrosine Kinases

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.4276 mL	12.138 mL	24.276 mL
5 mM	0.4855 mL	2.4276 mL	4.8552 mL
10 mM	0.2428 mL	1.2138 mL	2.4276 mL
50 mM	0.0486 mL	0.2428 mL	0.4855 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481