


c-Myc ligand 1

Chemical Properties

CAS No. :	1942834-87-8
Formula:	C13H10N2O2
Molecular Weight:	226.23
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	c-Mycligand 1 is a c-Myc inhibitor and functions as a ligand for target proteins in PROTAC (Proteolysis Targeting Chimeras). It is utilized in the synthesis of PROTAC c-Myc inhibitor 7.
Targets(IC50)	c-Myc,Ligands for Target Protein for PROTAC

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4203 mL	22.1014 mL	44.2028 mL
5 mM	0.8841 mL	4.4203 mL	8.8406 mL
10 mM	0.442 mL	2.2101 mL	4.4203 mL
50 mM	0.0884 mL	0.442 mL	0.8841 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481