

cRGD-PEG-Cy3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Powder: -20°C for 3 years | In solvent: -80°C for 1 year**  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	cRGD-PEG-Cy3 is a cRGD-PEG conjugate labeled with Cy3. The Cy3 fluorophore is commonly utilized in applications such as immunolabeling, nucleic acid labeling, fluorescence microscopy, and flow cytometry, with a peak emission wavelength of approximately 562-570 nm. cRGD can be used for lipid modification to enhance stability.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481