


1,2-Dihexanoyl-sn-glycero-3-PS sodium

Chemical Properties

CAS No. :	321862-85-5
Formula:	C ₁₈ H ₃₃ NNaO ₁₀ P
Molecular Weight:	477.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1,2-Dihexanoyl-sn-glycero-3-PS sodium is a phospholipid featuring hexanoic acid chains at both the sn-1 and sn-2 positions.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0946 mL	10.473 mL	20.9459 mL
5 mM	0.4189 mL	2.0946 mL	4.1892 mL
10 mM	0.2095 mL	1.0473 mL	2.0946 mL
50 mM	0.0419 mL	0.2095 mL	0.4189 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481