

## Eloralintide

## Chemical Properties

CAS No. :	2883634-40-8
Formula:	C201H319N49O65S2
Molecular Weight:	4526.1
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>

## Biological Description

Description	Eloralintide (LY 3841136) is an AMYR agonist, anticipated for research in type 2 diabetes and obesity.
Targets(IC50)	Amylin Receptor

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.2209 mL	1.1047 mL	2.2094 mL
5 mM	0.0442 mL	0.2209 mL	0.4419 mL
10 mM	0.0221 mL	0.1105 mL	0.2209 mL
50 mM	0.0044 mL	0.0221 mL	0.0442 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481