


Antitumor agent-160

Chemical Properties

CAS No. :	1262489-36-0
Formula:	C ₂₁ H ₃₀ O ₃
Molecular Weight:	330.46
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Antitumor agent-160 (compound 1a) is a derivative of the antifungal agent PHENALENONE and exhibits antitumor activity.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0261 mL	15.1304 mL	30.2608 mL
5 mM	0.6052 mL	3.0261 mL	6.0522 mL
10 mM	0.3026 mL	1.513 mL	3.0261 mL
50 mM	0.0605 mL	0.3026 mL	0.6052 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481