


1-Myristoyl-2-Linoleoyl-sn-glycero-3-PC

Chemical Properties

CAS No. :	92345-33-0
Formula:	C40H76NO8P
Molecular Weight:	730.01
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	1-Myristoyl-2-Linoleoyl-sn-glycero-3-PC is a phospholipid identified in human plasma, characterized by the presence of Myristic acid and Linoleic acid at the sn-1 and sn-2 positions, respectively.
Targets(IC50)	Phospholipase

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.3698 mL	6.8492 mL	13.6984 mL
5 mM	0.274 mL	1.3698 mL	2.7397 mL
10 mM	0.137 mL	0.6849 mL	1.3698 mL
50 mM	0.0274 mL	0.137 mL	0.274 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481